



Orthopedic Products

The Combination of comfort, technology, desing and functionality

Index



BodyCy Sports Knitted Supports	4-12
ThermoCy Neoprene Supports	13-23
OrthoCy Medical Supports	24-46
VariCy Varicose Socks	47-50
OneCy One Size Supports	51-56
KiddoCy Kids Supports	57-63
LifeCy Life Supports	64-66

ABOUT US

Through 17 years manufacturing experience in the soft orthopedic goods market in Turkey, we are providing our high-quality product line to more than 40 countries all around the world. We believe that the most valuable thing in life is health and we continue to work through the same consciousnesses from the first to the last step in our design and manufacturing process. Our goal is to grow sustainably in the global medical market with high-quality & innovative products and high-level customer satisfaction by putting the healthy and functional solutions on top of everything.

Morsa Orthopedics, the forerunner of the firsts in its field, presents a close to 200 product range to consumers under the name of MorsaCyberg brand, which is one of the first brands that come to mind in the medical orthopedics sector and has gained a serious place with quality and customer satisfaction.



Morsa Orthopedics knows that the most important thing in life is human health and continues to work in this direction with the same consciousness.

We are responsible for the environment as much as for our customers. Being the direct manufacturer helps us to control every single step of the sustainable manufacturing process including the eco-friendly material selection and to success the global challenge of going green in manufacturing.

Quality Management



QUALITY STANDARDS

All Products EU Medical Device Directive 93/42/EEC Annex I is according to the basic requirements.

QUALITY CERTIFICATES

ISO-13485:2016 ISO-9001:2015
CE Declaration of Conformity





BodyCy

Knitted Sports Orthopedic Products

The Combination of comfort,
technology, desing and functionality

KNITTED FABRICS

MorsaCyberg knitted braces feature three-dimensional knitting for a breathable, anti-allergenic, comfortable fit. As a supportive material, contoured with silicone pads to provide a high level of protection as well as compression, and stabilization during sports activities.

Knitted Fabrics; highly flexible, with breathable fabrics and compressive properties.
The MorsaCyberg Knitting products is among the high-tech machines are manufactured in technical textiles.

Ref-90.802



BodyCy Knitted Shoulder Support

- Provides support and warmth to the shoulder joints.
- The silicone pads provide full protection to the external effects.
- Made of highly elastic air permeable and compressive knitted fabric.

	S	34-38 cm	COLOR		Nude
	M	38-44 cm			Grey
	L	44-50 cm			Black
	XL	50-58 cm			
	XXL	58-66 cm			

BodyCy Knitted Elbow Support

Ref-50.304

- Used in the treatment of the tennis elbow syndrome and epicondylitis.
- The silicone pads provide support to flexor tendons.
- Made of highly elastic, air permeable and compressive knitted fabric.



	S	22-25 cm	COLOR		Nude
	M	25-28 cm			Grey
	L	28-31 cm			Black
	XL	31-34 cm			
	XXL	34-37 cm			

BodyCy Knitted Warm Elbow Support

Ref-180.102

- Provides warmth and support to the elbow joints.
- Soft knitted wickable and quick drying material.



	S	22-25 cm	COLOR		White
	M	25-28 cm			Blue
	L	28-31 cm			Claret Red
	XL	31-34 cm			
	XXL	34-37 cm			

Ref-60.419



BodyCy Knitted Wrist Support

- Stabilizes the wrist with the splint and provides immobilization.
- Used in the treatment process of carpal tunnel syndrome and wrist pain.
- Adjustable hook and loop tape provides an easy application.
- Made of highly elastic air permeable and compressive knitted fabric.

	S	15-17 cm	Color		Nude
	M	17-19 cm			Grey
	L	19-21 cm			Black
	XL	21-23 cm			
	XXL	23-25 cm			



Ref-60.421



BodyCy Knitted Wrist & Thumb Support

- Stabilizes the carpometacarpal and metacarpophalangeal joints of the thumb.
- Enables immobilization of the wrist in conditions such as carpal tunnel syndrome and painful wrist.
- Adjustable hook and loop tape provides an easy application.
- Made of highly elastic, air permeable and compressive knitted fabric.

	S	15-17 cm	Color		Nude
	M	17-19 cm			Grey
	L	19-21 cm			Black
	XL	21-23 cm			
	XXL	23-25 cm			



Ref-60.422



BodyCy Knitted Wrist Support

- Stabilizes the wrist with the splint and provides immobilization.
- Used in the treatment process of carpal tunnel syndrome and wrist pain.
- Made of highly elastic air permeable and compressive knitted fabric.

	S	15-17 cm	Color		Nude
	M	17-19 cm			Grey
	L	19-21 cm			Black
	XL	21-23 cm			
	XXL	23-25 cm			

Ref-180.100



BodyCy Knitted Warm Wrist Support

- The wrist is effective in situations that require heat and temperature.
- Soft knitted wickable and quick drying material.

C O L O R		White
		Blue
		Claret Red

	S	15-17 cm
	M	17-19 cm
	L	19-21 cm
	XL	21-23 cm
	XXL	23-25 cm

Ref-180.104



BodyCy Knitted Warm Corset

- Provides warmth and support to the back and abdominal area.
- Soft knitted wickable and quick drying material.

C O L O R		White
		Blue
		Claret Red

	S	70-80 cm
	M	80-90 cm
	L	90-100 cm
	XL	100-110 cm
	XXL	110-120 cm

Ref-40.402



BodyCy Knitted Lumbosacral Corset with Silicone Pad

- The silicone pad provides support to the lower back (lumbago) area.
- Used in the treatment process of facet joint disorders, sports injury, low back sprain, lumbar herniated disc, sacroiliac joint dysfunction, spondylolisthesis, osteoarthritis, back muscle strain, ligament sprain.
- Made of highly elastic air permeable and compressive knitted fabric



Color		Nude
		Grey
		Black


	S	70-80 cm
	M	80-90 cm
	L	90-100 cm
	XL	100-110 cm
	XXL	110-120 cm

BodyCy Knitted Lumbosacral Corset with Extra Belt

- Used in the treatment process of facet joint disorders, sports injury, low back sprain, lumbar herniated disc, sacroiliac joint dysfunction, back muscle strain and ligament sprain.
- The ergonomically shaped plastic plate provides high level of support to the lower back (lumbago) area.
- Made of highly elastic, air permeable and compressive knitted fabric.



Ref-40.404



Color		Nude
		Grey
		Black

	S	70-80 cm
	M	80-90 cm
	L	90-100 cm
	XL	100-110 cm
	XXL	110-120 cm

Ref-20.401



BodyCy Knitted Knee Support

- Used in the treatment process of osteoarthritis and meniscus lesions.
- Donut-Shaped silicone pad protects and stabilizes the knee.
- Made of highly elastic, air permeable and compressive knitted fabric.

S	34-36 cm	C O L O R		Nude
M	36-38 cm			Grey
L	38-41 cm			Black
XL	41-44 cm			
XXL	44-47 cm			

Ref-20.402



BodyCy Knitted Knee Support With Flexible Stays

- Used in the treatment of femoral pains, osteoarthritis and meniscus lesions.
- Donut Shaped silicone pad protects and stabilizes the knee.
- Flexible bonings support the medial and lateral ligaments of the knee.
- Made of highly elastic, air permeable and compressive knitted fabric.

S	34-36 cm	C O L O R		Nude
M	36-38 cm			Grey
L	38-41 cm			Black
XL	41-44 cm			
XXL	44-47 cm			

Ref-20.405



BodyCy Knitted Knee Support with Hinges

- Hinges on the both sides provide support to the knee in case of collateral traumas.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- Donut-Shaped silicone pad protects and stabilizes the knee.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Made of highly elastic, air permeable and compressive knitted fabric.





S	34-36 cm	C O L O R		Nude
M	36-38 cm			Grey
L	38-41 cm			Black
XL	41-44 cm			
XXL	44-47 cm			

Ref-20.403



BodyCy Knitted Knee Support with Hinges (Half Open)

- Provides support to the knee joint in case of strain of the knee ligaments.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- The steel springs on both of the side and perfect support to the lateral ligaments.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Half open front allows easy put on and take off.
- Made of highly elastic, air permeable and compressive knitted fabric.





	S	34-36 cm	C o l o r		Nude
	M	36-38 cm			Grey
	L	38-41 cm			Black
	XL	41-44 cm			
	XXL	44-47 cm			

Ref-20.404



BodyCy Knitted Knee Support with Hinges (Full Open)

- Provides support to the knee joint in case of strain of the knee ligaments.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- The steel springs on the sides provide perfect support to the lateral ligaments.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Full open front allows easy put on and take off.
- Made of highly elastic, air permeable and compressive knitted fabric.





	S	34-36 cm	C o l o r		Nude
	M	36-38 cm			Grey
	L	38-41 cm			Black
	XL	41-44 cm			
	XXL	44-47 cm			

BodyCy Knitted Warm Knee Support

- Provides warmth and support to the knee.
- Soft knitted wickable and quick drying material.

Ref-180.101



	S	34-36 cm	C o l o r		White
	M	36-38 cm			Blue
	L	38-41 cm			Claret Red
	XL	41-44 cm			
	XXL	44-47 cm			

Ref-90.441



BodyCy Knitted Ankle Support with Silicone Pad

- The silicone pads provide support to ankle.
- Used in the treatment process of edema.
- Increases blood circulation during strains.
- Made of highly elastic air permeable and compressive knitted fabric.

S	24-26 cm	C o l o r		Nude
M	26-28 cm			Grey
L	28-31 cm			Black
XL	31-34 cm			
XXL	34-37 cm			

Ref-90.442



Right / Left

BodyCy Knitted Achilles Ankle Support

- Provides support with silicone pads.
- Relieves the pain on achilles tendon ruptures and tendinitis edema.
- The pad placed over the Achilles tendon reduce the stress and the pad for the heel reduces the tension on the tendon.
- Made of highly elastic, air permeable and compressive knitted fabric.

S	24-26 cm	C o l o r		Nude
M	26-28 cm			Grey
L	28-31 cm			Black
XL	31-34 cm			
XXL	34-37 cm			

Ref-180.103




BodyCy Knitted Warm Ankle Support

- Provides warmth and support to the ankle ligament.
- It made from fabric with special woven acrylic and elastic fabric.

S	24-26 cm	C o l o r		White
M	26-28 cm			Blue
L	28-31 cm			Claret Red
XL	31-34 cm			
XXL	34-37 cm			

Ref-30.801



BodyCy Wool Shoulder Support

- Provides warmth to the shoulder joint.
- Provides light support in sport activities and in chronic pain conditions.
- Made of acrylic fabric.

	S	34-38 cm	C o l o r	White
	M	38-44 cm		
	L	44-50 cm		
	XL	50-58 cm		
	XXL	58-66 cm		

Ref-40.001



BodyCy Wool Corset

- Provides warmth to the abdominal and lumbar area by maintaining body temperature.
- Made of acrylic wool fabric.

	S	70-80 cm	C o l o r	White
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-20.001



BodyCy Wool Knee Support

- Provides warm and support to the knee.
- Made of acrylic wool fabric.

Ref-40.002



BodyCy Wool Corset (Open Front)

- Provides warmth to the abdominal and lumbar area by maintaining body temperature.
- Adjustable with front of velcro cover.
- Made of acrylic wool fabric.

	S	70-80 cm	C o l o r	White
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

	S	34-36 cm	C o l o r	White
	M	36-38 cm		
	L	38-41 cm		
	XL	41-44 cm		
	XXL	44-47 cm		



ThermoCy

Neopren

Orthopedic Products

The Combination of comfort,
technology, desing and functionality

NEOPRENE FABRICS


Neoprene; synthetic rubber is a product of the family. Is obtained by polymerization of one kind of rubber. Neoprene is widely used in underwater diver's clothes. MorsaCyberg neoprene knee supports and some other products are perferred in the perforated fabrics. (The perforated fabrics is prevents excessive sweating) we use neoprene fabrics in our products, plush coating on the inside of the anti-allergic and provides optimum comfort, the outer surface of the jersey fabric, provides maximum durability and flexibility.

MorsaCyberg Neoprene braces are made of breathable neoprene fabric with the channels on it which allow body heat and sweat go through as water vapor.

ThermoCy Knee Support Closed Patella

- Provides Support and warmth to the knee.
- Used in the treatment process of soft tissue injuries.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-20.100



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

C
o
l
o
r

Gray


Patterned Gray

Ref-20.101



ThermoCy Knee Support Open Patella

- Provides Support and warmth to the knee.
- Used in the treatment process of osteoarthritis and meniscus lesions.
- Donut Shaped silicone pad process and stabilizes the knee.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

C
o
l
o
r


Gray

Patterned Gray

ThermoCy Knee Support with Flexible Stays (Short)

- Provides Support and warmth to the knee.
- Used in the treatment process of osteoarthritis, ligament injuries, other knee joint traumas and meniscus lesions.
- Donut Shaped silicone pad process and stabilizes the knee.
- Flexible steel springs on the both sides provide convenient support to the knee in case of collateral traumas.
- Made of dual laminated perforated neoprene which allows air through the channels.

Ref-20.102



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

C
o
l
o
r

Gray


Patterned Gray

Ref-20.103




ThermoCy Knee Support With Flexible Stays (Long)

- Provides Support and warmth to the knee.
- Used in the treatment process of osteoarthritis , ligament injuries, other knee joint traumas and meniscus lesions.
- Donut Shaped silicone pad process and stabilizes the knee.
- Flexible steel springs on the both sides provide convenient support to the knee in case of collateral traumas.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

Color




Ref-20.104




ThermoCy Knee Support With Flexible Stays (Half Open)

- Provides strong support to the knee for arthritis, osteoarthritis, tendonitis,ligament ,sprains and strains.
- Support for Post-op or trauma for protection and rehabilitation purposes.
- The steel springs on either side designed to provide maximum support to the lateral ligaments.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Half open front part allows easy put on and take off.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

Color




Ref-20.015



ThermoCy Knee Support With Hinges

- Provides support to the knee joint in case of strain of the knee ligaments.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- Donut Shaped silicone pad process and stabilizes the knee.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Made of dual laminated perforated neoprene which allows air through the channels.
- Made of air - porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm


Color



ThermoCy Knee Support with Hinges

Ref-20.106

- Provides support to the knee joint in case of strain of the knee ligaments.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- Donut Shaped silicone pad process and stabilizes the knee.
- The Hinged on either side designed to provide maximum support to the lateral ligaments.
- With the upper and lower hook and loop straps, it is aimed to prevent slipping down.
- Half open front part allows easy put on and take off.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm


Color		Gray
		Patterned Gray



ThermoCy Knee Support with Hinges (Full Open)

Ref-20.107

- Provides support to the knee joint in case of strain of the knee ligaments.
- Used in the ligament injuries process of soft tissue lesions and other knee joint treatments.
- Donut Shaped silicone pad process and stabilizes the knee.
- Hinges on either side designed to provide maximum support to the knee
- Upper and lower hook and loop bandages prevent slipping down
- Full open front part allows easy put on and take off.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm


Color		Gray
		Patterned Gray



ThermoCy Knee Support with Patella Pad

Ref-20.108

- Provides protection for patella and the knee against impact especially during sports activities.
- The product functions well with patella and soft tissue lesions around the patella.
- The soft pad covering the patella absorbs shocks to the knee.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

Color		Gray
		Patterned Gray



ThermoCy Knee Support with Cross Ligament

Ref-20.109

- Provides Support to the anterior cruciate ligament with its cross straps.
- The cross-velcro strap, which can be adjusted in both directions, prevents excessive stretching of the cross-links prevents not slipping back and forth.
- The steel springs on either side designed to provide maximum support to the lateral ligaments.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm	Color	<input type="checkbox"/> Gray <input checked="" type="checkbox"/> Patterned Gray
M	36-38 cm		
L	38-41 cm		
XL	41-44 cm		
XXL	44-47 cm		

Ref-20.111



ThermoCy Knee Tendon Support

- Used in treatment process of Osgood-Schlatter's syndrome and patellar tendinitis.
- Thanks to the air channels in the neoprene, it is aimed to reduce excessive perspiration without losing body temperature.

S	34-36 cm	Color	<input type="checkbox"/> Gray <input checked="" type="checkbox"/> Patterned Gray
M	36-38 cm		
L	38-41 cm		
XL	41-44 cm		
XXL	44-47 cm		

ThermoCy Knee Support with Controlled Mechanism

Ref-83.013

- Used in the treatment process of soft tissue and ligament injuries and other knee joint traumas that does not need ACL and PCL operations.
- Inhibits the movement to desired levels when knee control and protection is required owing to the adjustable angle bilateral joint.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.



S	34-36 cm	Color	<input type="checkbox"/> Gray <input checked="" type="checkbox"/> Patterned Gray
M	36-38 cm		
L	38-41 cm		
XL	41-44 cm		
XXL	44-47 cm		

Ref-40.201



ThermoCy Corset with Stays

- Provides warmth and support to the lumbar area.
- Used in the treatment process of lumbar and rheumatismal pain.
- Appropriate to use during exercise and training.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	70-80 cm	C o l o r	Gray
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		
			Patterned Gray	

ThermoCy Corset With Belt

- Provides warmth and support to the lumbar area.
- Used in the treatment process of lumbar and rheumatismal pain
- Appropriate to use during exercise and training.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration

Ref-40.202




	S	70-80 cm	C o l o r	Gray
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		
			Patterned Gray	

ThermoCy Abdominal Corset

- Provides warmth and support to the abdominal area.
- Appropriate to use during exercise and training.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-40.203



	S	70-80 cm	C o l o r	Gray
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		
			Patterned Gray	

Ref-50.301



ThermoCy Elbow Support

- Provides warmth and support to the elbow joint.
- Used in the treatment process of tennis elbow syndrome and epichondylitis.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	22-25 cm	C O L O R		Gray
	M	25-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			

ThermoCy Elbow Support With Pad

- Provides support to the elbow area and against impact.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-50.302



	S	22-25 cm	C O L O R		Gray
	M	25-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			

ThermoCy Wrist Support

- Provides support by restricting the movement of the hand wrist for Carpal Tunnel Syndrome, wrist fractures, wrist sprains and wrist tendonitis.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-60.601



	S	15-17 cm	C O L O R		Gray
	M	17-19 cm			Patterned Gray
	L	19-21 cm			
	XL	21-23 cm			
	XXL	23-25 cm			

Ref-60.604



ThermoCy Wrist & Thumb Support

- Provides excellent immobilization to the thump joint and gives an external support.
- Functions well with De Quarvain syndrome and thumb ligament problems.
- Stabilizes the wrist thanks to the aluminum thumb support.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	15-17 cm	C O L O R		Gray
	M	17-19 cm			Patterned Gray
	L	19-21 cm			
	XL	21-23 cm			
	XXL	23-25 cm			

ThermoCy Thigh Support

- Provides femur protection and also support quadriceps and adductor muscle groups with head and compression.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-70.500



	S	48-52 cm	C O L O R		Gray
	M	52-57 cm			Patterned Gray
	L	57-62 cm			
	XL	62-70 cm			
	XXL	70-78 cm			

ThermoCy Calf Support

- Provides protection, compression and warmth to the tibial area in case of muscular leisons.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-70.501



	S	34-36 cm	C O L O R		Gray
	M	36-38 cm			Patterned Gray
	L	38-41 cm			
	XL	41-44 cm			
	XXL	44-47 cm			

Right / Left

Ref-70.503



ThermoCy Hip Support


- Provides protection warmth and support to the soft tissues around the hip joint.
- Supports the hip joint with a plastic or metal hinge placed in a pocket located on the hip joint.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Right / Left	Color	Sizing		
			Gray	
		Patterned Gray		
	S	70-80 cm		48-52 cm
	M	80-90 cm		52-57 cm
	L	90-100 cm		57-62 cm
	XL	100-110 cm		62-70 cm
	XXL	110-120 cm		70-78 cm

ThermoCy Shoulder Support

- Provides Protection and warmth to the shoulder joints.
- Used in the treatment of chronic pains that are caused by stress or sport activities.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-30.800



Color	Sizing	
		Gray
	Patterned Gray	
	S	34-38 cm
	M	38-44 cm
	L	44-50 cm
	XL	50-58 cm
	XXL	58-66 cm



ThermoCy Sport Swimsuit

- Provides ease of use with zip on front fasteners and thanks to the elastic fabric used under the chest.

Ref-70.702



	S	70-80 cm	C O L O R	
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

ThermoCy Sport Tights (Short)

- Provides protection to the exercise injuries.
- Helps to lose weight from the waists, hip, and abdominal areas.
- Tighten up and smoothen the skin.
- Helps to get rid of toxins and cellulite.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-70.700



Ref-70.700



Ref-70.700



	S	70-80 cm	C O L O R	
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

ThermoCy Sport Tights (Long)

- Provides protection to the exercise injuries.
- Helps to lose weight from the waists, hip, and abdominal areas.
- Tighten up and smoothen the skin.
- Helps to get rid of toxins and cellulite.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

Ref-70.701



Ref-70.701



Ref-70.701



	S	70-80 cm	C O L O R	
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-90.401



ThermoCy Ankle Support

- Provides warmth and support to the stabilization for the mild sprains.
- Provides stronger support with the cross strap.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	24-26 cm	Color		Gray
	M	26-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			

Ref-90.403



ThermoCy Ankle Support With Pad

- Provides warmth and support to the stabilization for the mild sprains.
- Provides massage effect with pad to the ankle joint and helps to remove hematoma and edema.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	24-26 cm	Color		Gray
	M	26-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			

Ref-90.404



ThermoCy Ankle Support With Reinforcement (Laced)

- Provides support and warmth to the ankle in cases of mild sprains.
- Effective with habitual ankle sprains, chronic ankle pains and can be used as a protector with contact sports.
- Thanks to two specially designed plastic plates located on the malleols, it provides strong support to the ankle.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	24-26 cm	Color		Gray
	M	26-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			

Ref-90.408



ThermoCy Ankle Support With Reinforcement (Laced)

- Provides stronger stabilization to the ankle sprains.
- Effective with habitual ankle sprains, chronic ankle pains and can be used as a protector with contact sports.
- Thanks to two specially designed plastic plates located on the malleols, it provides strong support to the ankle.
- Made of air-porous neoprene fabric with aimed to reduce excessive perspiration.

	S	24-26 cm	Color		Gray
	M	26-28 cm			Patterned Gray
	L	28-31 cm			
	XL	31-34 cm			
	XXL	34-37 cm			



Morsa
Cyberg

OrthoCy

Medical Orthopedic

The Combination of
comfort, technology,
desing and functionality



NeckCy Cervical Orthosis

- Used in the treatment process of injuries on the cervical spine by supporting the chin.
- Constrains movements of the neck following neck traumas and cervical disc hernia.
- Neck height is adjustable.
- Made of PVC material.

	S	34-36 cm	Color	White
	M	36-38 cm		
	L	38-40 cm		
	XL	40-42 cm		
	XXL	42-44 cm		

NeckCy Cervical Miorva Orthosis (CTO)

- Provides traction to the neck spine by immobilizing the head and the trunk.
- Used in the treatment process of serious neck traumas and post neck operations.
- Enables the spine to be secured and fixed in the correct position.
- Hook and Loop bandages enable easy adjustment.
- The frame is made of polyethylene material.



Ref-11.105



	White		
	S	A 34-36 cm	B 70-90 cm
	M	A 36-38 cm	B 90-110 cm
	L	A 38-40 cm	B 110-130 cm
	XL	-	-
	XXL	-	-

Ref-11.106



NeckCy Cervical Orthosis (First Aid)

- First Aid necessary to avoid the compression of the cervical for a polytraumatized patient during moving or transportation.
- Designed for use in first aid.
- Made of PVC material.


Adult - Child		Adjustable	Color	White

Ref-11.107



NeckCy Cervical Orthosis (Traction)

- Provides excellent stabilisation following neck traumas and cervical disc hernia.
- Thanks to its special design, it offers optimum comfort.
- It is made of polyethylene material covered with sponge fabric laminated sponge.




S	34-36 cm	Color	<input type="checkbox"/> Nude
M	36-38cm		
L	38-40cm		
XL	-		
XXL	-		

Ref-11.109



NeckCy Cervical Orthosis (Wire Frame Traction)

- Provides excellent stabilization following neck traumas and cervical disc hernia.
- Thanks to its special design, it offers optimum comfort.
- The front part is made of aluminum frame and back part is made of polyethylene material.




S	34-36 cm	Color	<input type="checkbox"/> Gray
M	36-38cm		
L	38-40 cm		
XL	-		
XXL	-		

Ref-12.102



NeckCy Cervical Orthosis (Foam)

- Provides limited cervical motion control and relieves the mild neck pain by limiting the cervical spine motions.
- Used in the treatment process of the disk hernia cervical osteoarthritis and cervical muscle pains.
- Hook and pile closure in back for quick and easy fit.
- Made of plastazote (compressed foam) material in thickness.



S	34-36 cm	Color	<input type="checkbox"/> White <input checked="" type="checkbox"/> Darkblue
M	36-38 cm		
L	38-40 cm		
XL	40-42 cm		
XXL	42-44 cm		

Ref-12.103



NeckCy Cervical Orthosis (Sponge) 7.5 cm - 8.5 cm - 10 cm

- Provides limits cervical spine motion to varying degrees.
- Used in the treatment of a disc hernia cervical osteoarthritis and cervical muscle pains.
- Hook and pile closure in back for quick and easy fit.
- Made of density sponge in thickness.



	S	34-36 cm	C O L O R		Grey
	M	36-38 cm			Yellow
	L	38-40 cm			Mustard Yellow
	XL	40-42 cm			Pink
	XXL	42-44 cm			Green
					White
				Darkblue	

Ref-12.104



NeckCy Cervical Orthosis (Philadelphia Type)

- Provides support and helps to immobilise the neck following soft tissue or disc trauma.
- Used to stabilize minor cervical spinal fractures and reduce pain associated with a muscle strain on the back.
- Supports a weak or injured neck, and helps with pain and instability.
- Lightweight and flexible material for comfort and durability.
- Suitable for X-RAY and MR monitoring.
- Made of plastozote (compressed foam) material in thickness.

	S	34-36 cm	C O L O R		Nude
	M	36-38cm			
	L	38-40cm			
	XL	-			
	XXL	-			

NeckCy Cervical Orthosis (Trachostomy - Open)

- Provides support and helps to immobilise the neck following soft tissue or disc trauma.
- Used to stabilize minor cervical spinal fractures and reduce pain associated with a muscle strain on the back.
- Supports a weak or injured neck, and helps with pain and instability.
- Light weight and flexible material for comfort and durability.
- Suitable for X-RAY and MR monitoring.
- Made of plastozote (compressed foam) material in thickness.

Ref-12.108



	S	34-36 cm	C O L O R		Nude
	M	36-38 cm			
	L	38-40 cm			
	XL	-			
	XXL	-			

Ref-13.110



ProCy Nose Protective Mask

- Provides protection pre/post nose operation or during sports activities.
- Protects from further nose injury.
- Helps protect user from minor bumps.
- Special pad supports on the forehead and cheek area.
- No vision obstruction.
- Easily adjusted by elastic bands.
- Made of acrylic material.



S	22 / 22,5 cm	Color
M	23 / 23,5 cm	
L	24 / 24,5 cm	
XL	-	
XXL	-	

SoftCy Arm Sling (Lux)

- Provides support for broken wrist, upper arm and fore arm injuries.
- Stabilizes the arm in the desired position to immobilize the arm-shoulder.
- Comes with an adjustable neck-shoulder strap and a waist belt.
- Made of sponge laminated with terry cloth in the inside and flocculent fabric in the outside.

Ref-30.301



	S	34-38 cm	Color	Grey Darkblue
	M	38-44 cm		
	L	44-50 cm		
	XL	50-58 cm		
	XXL	58-66 cm		

Ref-30.302



SoftCy Arm Sling (Mesh fabric)

- Provides support for broken wrist, upper arm and forearm injuries.
- Stabilizes the arm in the desired position to immobilize the arm-shoulder.
- Made of sponge fabric.

	S	34-38 cm	Color	Grey Darkblue
	M	38-44 cm		
	L	44-50 cm		
	XL	50-58 cm		
	XXL	58-66 cm		

Ref-30.303



SoftCy Arm Sling (Gabardin Fabric)

- Provides support for broken wrist, upper arm and forearm injuries.
- Stabilizes the arm in the desired position to immobilize the arm-shoulder.
- Adjustable neck and waist strap.
- Made of gabardin fabric.

	S	34-38 cm	C O L O R	<input type="checkbox"/> White
	M	38-44 cm		<input type="checkbox"/> Dark Blue
	L	44-50 cm		
	XL	50-58 cm		
	XXL	58-66 cm		

Ref-30.304



SoftCy Shoulder Immobilizer

- Provides excellent immobilization to the shoulder and upper arm in the treatment process of humerus fractures and shoulder dislocations.
- Made of sponge with laminated fabric.

	S	70-80 cm	C O L O R	<input type="checkbox"/> Grey
	M	80-90 cm		<input type="checkbox"/> Dark Blue
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

SoftCy Arm Sling With Cushion 30 Deg.

- Provides immobilization limits the motion of the arm at 30 degrees of abduction.
- Used in the treatment of the post rotator cuff, habitual dislocations and subluxation operation.
- The design of the cushion allows easy use of either arm.
- Made of laminated fabric and sponge cushion.

Ref-30.305



	Adjustable size	C O L O R	<input type="checkbox"/> Grey
			<input type="checkbox"/> Dark Blue

Ref-30.306



OrthoCy Arm Sling With Cushion 45-60 Deg.

- Provides immobilization limits the motion of the arm at 45-60 degrees of abduction.
- Used in the treatment of the post rotator cuff, habitual dislocations and subluxation operation.
- The design of the cushion allows easy use of both arm.
- Made of laminated fabric and sponge cushion.

Adjustable size

COLOR		Grey
		Dark Blue

SoftCy Clavicle Bandage


- It can replace a plaster to stabilize the clavicle bone in cases of clavicle fractures.
- The one-piece adjustable bandage is designed- for the treatment of clavicle fractures.
- Made of laminated sponge in thickness.

Ref-30.311



COLOR	S	70-80 cm		Grey
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		
				Dark Blue

Ref-31.800



OrthoCy Humerus Brace

- Humerus brace used for stabilization and protection in upper extremity fractures surgery and post-plaster shoulder fractures.
- Made of thermoplastic material.



Right- Left

COLOR	S	32-35 cm		Grey
	M	36-39 cm		
	L	40-43 cm		
	XL	-		
	XXL	-		
				Dark Blue

LumbaCy Corset Posturex

- Provides support in case of postural deficiency and mild kyphosis.
- Elastic shoulder bandages enable upright position and help patients upright sitting habits.
- Made of gutted elastic fabric that is woven with cotton yarn.

Ref-40.161



	S	70-80 cm	Color	Nude
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-40.162



LumbaCy Corset Posturex For Pregnancy

- Provides support to the waist and back muscles at the same time and the abdomen used to prevent sagging.
- Used to corrective effect in prenatal period, postural disorders and mild degrees of kyphosis.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S	70-80 cm	Color	Nude
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-40.250



LumbaCy Corset Lumbosacral Without Belt

- Provides support to the lower back.
- Used in the treatment process of sacroiliac and lumbosacral pain.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Back height is 26 cm.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S	70-80 cm	Color	Nude
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-40.260



LumboCy Corsel Lumbosacral With Belt (26 cm)

- Provides support to the lower back.
- Used in the treatment process of sacroiliac and lumbo sacral pain.
- With additional bandages, stronger support to the lumbar muscles is acquired.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Back height is 26 cm.
- Made of gutted elastic fabric that is woven with cotton yarn.
- With additional bandages, stronger support to the lumbar muscles is acquired.

	S 70-80 cm	C O L O R	Nude
	M 80-90 cm		
	L 90-100 cm		
	XL 100-110 cm		
	XXL 110-120 cm		

LumboCy Corsel Lumbosacral (Adjustable)

- Provides support to the lower back.
- Used in the treatment process of sacroiliac and lumbo sacral pain.
- Thanks to all kinds of bandage can be adjusted in the body structure that provides better support.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Back height is 26 cm.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S 70-80 cm	C O L O R	Nude
	M 80-90 cm		
	L 90-100 cm		
	XL 100-110 cm		
	XXL 110-120 cm		

Ref-40.270



Ref-40.320



LumboCy Corsel Lumbosacral With Belt (32 cm)

- Provides support to the lower back.
- Used in the treatment process of sacroiliac and lumbo sacral pain.
- With additional bandages, stronger support to the lumbar muscles is acquired.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Back height is 32 cm.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S 70-80 cm	C O L O R	Nude
	M 80-90 cm		
	L 90-100 cm		
	XL 100-110 cm		
	XXL 110-120 cm		

Ref-40.330



LumbaCy Corset Thorax Lumber

- Provides support to correct dorsal spine and used in the treatment process of kphosis.
- Support to the stabilizes hyperextension in vertebral fractures of lumbar and thoracic spine.
- Additional straps provide stronger support to the waist and back muscles.
- Removable steel stays at the back enable easier adjustment according to the form of the vertebra.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S	70-80 cm	C o l o r	Nude
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

LumbaCy Corset Abdominal

- Provides support to the abdominal area following surgeries.
- Prevents sagging in the abdominal area following pregnancy.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Made of gutted elastic fabric that is woven with cotton yarn.

Ref-40.420



	S	70-80 cm	C o l o r	Nude
	M	80-90 cm		
	L	90-100 cm		
	XL	100-110 cm		
	XXL	110-120 cm		

Ref-40.421



LumbaCy Corset Sternum

- Provides stabilization in treatment of fractures of the sternum.
- Supports sternum after an open heart or thoracic surgery.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Made of gutted elastic fabric that is woven with cotton yarn.

	S	80-90 cm	C o l o r	Nude
	M	90-100 cm		
	L	100-110 cm		
	XL	110-120 cm		
	XXL	120-130 cm		

Ref.40.422



LumbaCy Corset Umbilical Hernia

- Abdominal hernias are supported by applying pressure through the special pad and prevents the progression of hernia.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Made of gutted elastic fabric that is woven with cotton yarn.


	S 70-80 cm	C O L O R	<input type="checkbox"/> Nude
	M 80-90 cm		
	L 90-100 cm		
	XL 100-110 cm		
	XXL 110-120 cm		

LumbaCy Corset For Pregnancy

- Provides support to the sacroiliac muscles to prevent lower back problems.
- Appropriate to use during and after pregnancy period to provide support to the sacroiliac area.
- Helps to prevent sagging on the abdominal area after pregnancy.
- Thanks to the weaving system of the fabric with air channels, perspiration is minimized.
- Made of gutted elastic fabric that is woven with cotton yarn.

Ref.40.439




	S 70-80 cm	C O L O R	<input type="checkbox"/> Nude
	M 80-90 cm		
	L 90-100 cm		
	XL 100-110 cm		
	XXL 110-120 cm		

Ref.41.400



OrthoCy Thoraco Lumbosacral Orthosis (TLSO)

- Its use to prevents lateral bending and rotation in the cases of hard trauma, in vertebral fractures and trunk flexion extension treatment period.
- Used in rehabilitation after surgery.
- Made of thermoplastic material.

	S 70-80 cm	C O L O R	<input type="checkbox"/> White
	M 80-90 cm		
	L 90-100 cm		
	XL -		
	XXL -		

Ref-44.550



OrthoCy Corset Jewel Hyperextension

- The hyperextension of vertebral column provides stabilization on the sagittal plane.
- Prevents lateral flexion with aluminum frame.
- Consists of a front part that is an aluminum frame and a pack part that is a plastic support laminated with sponge and fabric.



Adjustable size

COLOR

Blue

OrthoCy Corset Jewel Hyperextension

- The hyperextension of vertebral column provides stabilization on the sagittal plane.
- Prevent limits lateral flexion and rotation of the upper body.

Ref-44.555



S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	-
XXL	-

COLOR

White

Ref-53.015



OrthoCy Elbow Contracture Orthosis

- Enables Controlled mobilization of elbow joint.
- Used in the temporary immobilization following medial and lateral epicondyl fractures, elbow luxation / dislocation, elbow arthrolysis.



Adjustable size

COLOR

Grey
Black

OrthoCy Elbow Contracture Orthosis (Thermoplastic)

- Enables Controlled mobilization of elbow joint.
- Used in the temporary immobilization following medial and lateral epicondyl fractures, elbow luxation/dislocation, elbow arthrolysis.

Ref-53.016



Right- Left

	S	22-28 cm	Color	Grey
	M	28-34 cm		
	L	34-40 cm		
	XL	-		
	XXL	-		

Ref-60.015



OrthoCy Thumb Support

- Inhibits the thumb joint and provides external support.
- Used in De quarvain syndrome and thumb ligament problems.
- Made of specially designed elastic cotton fabric.
- Thanks to the aluminum rod inside, position of the thumb is adjustable.

	S	15-17 cm	Color	Nude Black
	M	17-19 cm		
	L	19-21 cm		
	XL	21-23 cm		
	XXL	23-25 cm		

Ref-60.016



OrthoCy Wrist Support (Short)

- Inhibits the thumb joint and provides external support.
- Used in De quarvain syndrome and thumb ligament problems.
- Made of specially designed elastic cotton fabric.
- Thanks to the aluminum rod inside, position of the thumb is adjustable.

Right - Left

	S	15-17 cm	Color	Nude Black
	M	17-19 cm		
	L	19-21 cm		
	XL	21-23 cm		
	XXL	23-25 cm		

Ref-60.019



OrthoCy Wrist Support (Long)

- Inhibits the thumb joint and provides external support.
- Used in De quarvain syndrome and thumb ligament problems.
- Made of specially designed elastic cotton fabric.
- Thanks to the aluminum rod inside, position of the thumb is adjustable.

Right - Left

	S	15-17 cm	Color	Nude Black
	M	17-19 cm		
	L	19-21 cm		
	XL	21-23 cm		
	XXL	23-25 cm		




Ref-60.021

OrthoCy Wrist & Thumb Support

- Stabilizes the carpometacarpal andmetacar pophalangeal joints of the thumb.
- Enables immobilization of the wrist in conditions such as carpal tunnel syndrome and painful wrist.
- Stabilizes the carpometacarpal andmetacar pophalangeal joints of the thumb.
- Enables immobilization of the wrist in conditions such as carpal tunnel syndrome and painful wrist.

Left - Right

	S	15-17 cm	Color	<input type="checkbox"/> Nude <input type="checkbox"/> Black
	M	17-19 cm		
	L	19-21 cm		
	XL	21-23 cm		
	XXL	23-25 cm		



Ref-61.001

OrthoCy Finger Support (Mallet Finger)

- This splint is used in the treatment of Mallet finger that Refers to the involuntary flexion of the distal phalanx of a finger caused by the disruption or tearing of its extensor tendon.
- Made of plastic material by making holes in order to provide better ventilation.

	S	4,5-5 cm	Color	<input type="checkbox"/> Nude
	M	5-6 cm		
	L	6-7 cm		
	XL	7-8 cm		
	XXL	8-9 cm		




Ref-61.031

OrthoCy Static Thumb Orthosis

- Stabilizes the carpometacarpal and meta carpophalangeal joints of the thumb.
- Made of specially designed elastic cotton.
- Made of polyethylene sheets in 3 mm thickness.
- Provide optimum comfort and cleanliness thanks to the removable holster.

Left - Right

	S	15-17 cm	Color	<input type="checkbox"/> Dark Blue <input type="checkbox"/> Grey
	M	17-19 cm		
	L	19-21 cm		
	XL	-		
	XXL	-		



Ref-61.033

OrthoCy Static Wrist Orthosis

- Enables immobilization of wrist.
- Prevents ulnar and radial rotations.
- Used in resting position following plaster treatment of fractures and painful positions of the wrist. Made of specially designed elastic cotton fabric.
- Made of polyethylene sheets in thickness.
- Provide optimum comfort and cleanliness thanks to the removable holster.

Left - Right

	S	15-17 cm	Color	<input type="checkbox"/> Dark Blue <input type="checkbox"/> Grey
	M	17-19 cm		
	L	19-21 cm		
	XL	-		
	XXL	-		



Ref-61.034

OrthoCy Static Hand Orthosis

- Provides a resting position for wrist by keeping wrist at 30 degrees' dorsiflexion and IP joints at mild flexion.
- Made of polyethylene sheets in 3 mm thickness.
- Provide optimum comfort and cleanliness thanks to the removable holster.

Left - Right

	S	15-17 cm	Color	<input type="checkbox"/> Dark Blue <input type="checkbox"/> Grey
	M	17-19 cm		
	L	19-21 cm		
	XL	-		
	XXL	-		



Ref-61.035

OrthoCy Static Hand-Finger Orthosis

- Provides a resting position for wrist by keeping wrist at 30 degrees' dorsiflexion and IP joints at mild flexion.
- Made of polyethylene sheets in 3 mm thickness.
- Provide optimum comfort and cleanliness thanks to the removable holster.

Left - Right

	S	15-17 cm	Color	<input type="checkbox"/> Dark Blue <input type="checkbox"/> Grey
	M	17-19 cm		
	L	19-21 cm		
	XL	-		
	XXL	-		



OrthoCy Dynamic Hand Orthosis

- Commonly used to exercise the fingers in cases of stiffness or inactivity following Cerebral Palsy, Distal Radius/ Ulna, Multiple Sclerosis, Muscular Dystrophy, Stroke and tendon & ligament operations.
- Easy to adjust
- Biomechanically correct
- Applicable with children and adults.

S	15-17 cm	Color	Dark Blue
M	17-19 cm		
L	19-21 cm		
XL	-		
XXL	-		
			Grey



OrthoCy Dynamic Hand Orthosis (Functional)

- Commonly used to exercise the fingers in cases of stiffness or inactivity following Cerebral Palsy, Distal Radius/Ulna, Multiple Sclerosis, Muscular Dystrophy, Stroke and tendon & ligament operations.

S	15-17 cm	Color	Dark Blue
M	17-19 cm		
L	19-21 cm		
XL	-		
XXL	-		
			Grey



OrthoCy Finger Support (Frog)

- Maintains both interphalangeal joints in natural functional position.
- Can be folded to accommodate various finger sizes.
- Made of foam laminated aluminum plate for optimum comfort.


S	4 - 6 cm	Color	Orange
M	6 - 8 cm		
L	8 - 10 cm		
XL	-		
XXL	-		



OrthoCy Finger Support (Baseball)

- Used to maintain interphalangeal joints in correct position following distal joint injuries.
- This splint is used to hold the finger in alignment
- Can be folded to accommodate various finger sizes.
- Made of foam laminated aluminum plate for optimum comfort.

S	4 - 6 cm	Color	Orange
M	6 - 8 cm		
L	8 - 10 cm		
XL	-		
XXL	-		



OrthoCy Finger Support (Absorbium Stick)

- Applicable as finger cot, baseball splint, finger protector, or for transverse fractures of the phalanges
- Made of foam laminated aluminum plate for optimum comfort.

30 cm		Color	Orange



OrthoCy Dynamic Finger Support

- Extension allows the finger immobilization.
- Finger fractures, tendon and nerve injuries, sprains hurt and used.

S	4 - 6 cm	Color	Nude
M	6 - 8 cm		
L	8 - 10 cm		
XL	-		
XXL	-		

Ref-70.101




Left - Right

OrthoCy Inguinal Rupture Belt

- Used for Inguinal hernia before and after surgery.
- Specially designed to provide pressure with silicone pad on the hernia.
- Easy to apply with adjustable waist and groin straps.

	S	Adjustable	C o l o r		Nude
	M				
	L				
	XL				
	XXL				

Ref-70.102



OrthoCy Inguinal Rupture Belt (Double)

- Used for Inguinal hernia before and after surgery.
- Specially designed to provide pressure with silicone pad on the hernia.
- Easy to apply with adjustable waist and groin straps.

	S	Adjustable	C o l o r		Nude
	M				
	L				
	XL				
	XXL				

OrthoCy Supporter (S) Suspensory

- In support of the testicles and is used to prevent sagging.

Ref-70.110



	S	70-80 cm	C o l o r		White
	M	80-90 cm			
	L	90-100 cm			
	XL	100-110 cm			
	XXL	110-120 cm			



OrthoCy Hip Abduction Orthosis with Functional Joint

- Post-operative stabilization after total hip replacement and revision.
- Prevents recurrent luxation.
- Abduction and Adduction angles can be adjusted.
- Flexion and extension angles can be adjusted in 20 degrees' increments (up to 120 degrees)
- Protects active transmission in the pelvic and inguinal area.


	A	B	Color
S	70-80 cm	48-52 cm	
M	80-90 cm	52-57 cm	Grey
L	90-100 cm	57-62 cm	
XL	-	-	
XXL	-	-	



OrthoCy Knee Immobilizer (45 cm) - (55 cm) - (65 cm) - (75 cm)

- Used in situations where full extension is needed in knee joint.
- Provides post-operative support following patellar, arthral and collateral ligament operations.
- Made with fabric laminated sponge in 6 mm thickness laminated with terycloth interiorly.
- Provides 0 extension to the knee with stays on both sides and at the back axis of the knee.

Ref-80.045	Adjustable	Color	Dark Blue
Ref-80.055			Grey
Ref-80.065			
Ref-80.075			



OrthoCy Knee Brace With Controlled Mechanism


- Inhibits the movement to desired levels when knee control and protection is require owing to the adjustable angle bilateral joint.
- Provides excellent support following ACL, PCL operations.
- Thanks to the adjustable bilateral hinges, it enables the desired flexion restriction.
- There are skin-friendly pads are removable and washable.



Adjustable	Color	Dark Blue
		Grey

OrthoCy Knee Brace with Controlled Mechanism (moving length)

Ref-83.012

- Inhibits the movement to desired levels when knee control and protection is require owing to the adjustable angle bilateral joint.
- Provides excellent support following ACL, PCL operations.
- Thanks to the adjustable bilateral hinges, it enables the desired flexion restriction.
- There are skin-friendly pads are removable and washable.




Adjustable	Color	 Dark Blue
		 Grey





OrthoCy Multi Knee Brace

Ref-83.014

- In cases of Gonarthrosis (Varus or Valgus deformities) helps to correct the knee axis where knee compartment has reduced.
- Provides excellent support in cases of ACL and PCL operations as well as other cases where control and protection is needed.
- Thanks to its specially designed aluminum frame, it maintains anatomical position of the knee.
- Thanks to the adjustable bilateral hinges, it enables the desired flexion restriction.



Color	S	34-36 cm	 Dark Blue
	M	36-38 cm	
	L	38-41 cm	
	XL	41-44 cm	
	XXL	44-47 cm	
			 Grey






Left - Right

OrthoCy Ankle Foot Orthosis (ReFlex AFU)

Ref-91.011

- In drop foot symptom, prevents friction by keeping the foot and ankle at the functional position
- Polypropylene is made of high-endurance and flexibility. Thanks to a specially designed mold can be easily applied into any kind of shoes.



Color	S	35-37 cm	 Dark Blue
	M	38-40 cm	
	L	40-44 cm	
	XL	-	
	XXL	-	
			 Grey



Left - Right



Ref-91.051

OrthoCy Ankle Foot Orthosis (Foot Guard)

- Enables resting of the foot and ankle in multi-positions.
- Prevents eschar wounds
- Made of polythene in 3 mm thickness.
- Provides comfort and cleanliness thanks to removable terrycloth pads.
- Thanks to the adjustable T shaped splint at the bottom, it enables adjustment of height and balance of the foot.

S	34-37 cm	Color	Darkblue
M	38-41 cm		Grey
L	42-45 cm		
XL	-		
XXL	-		



Ref-91.101

OrthoCy Ankle Support With Air Pad

- Functions as a stabilizer in cases of serious ankle sprains and chronic ankle pains and as a protector in contact sports.
- Also applicable in cases of ankle fractures and in the period after plastering process.
- Made with specially designed plastic moulds.
- Thanks to the air pads inside, it fits to the foot anatomically and provides maximum stabilisation.

Adjustable	Color	White
------------	-------	-------



Ref-91.102

OrthoCy Ankle Support With Plastic Hull

- With its rigid exterior, provides good stabilization in pronation and supination actions.
- Provides full support for ankle ligament injuries.
- The hard outer shell has been made with specially designed plastic moulds.
- Thanks to the soft pads inside, it fits to the foot anatomically and provides maximum stabilisation.

Adjustable	Color	Black
------------	-------	-------



Ref-93.001

OrthoCy Rom Walker

- Inhibits the movement to desired levels when ankle control and protection is required owing to the adjustable angle bilateral joint
- With the side angle adjustable ankle joints is provided in the control of the desired degree.

S	34-37 cm	Color	Black
M	38-41 cm		
L	42-45 cm		
XL	-		
XXL	-		



Ref-93.002

OrthoCy Achilles Shoes

- After conservative and surgical treatment of Achilles tendon rupture is used for protection or early functional treatment
- Thanks to each of the Achilles tendon in 10-degree heel risers may be desired position

S	34-37 cm	Color	Black
M	38-41 cm		
L	42-45 cm		
XL	-		
XXL	-		



Ref-96.001

OrthoCy Cast Shoes

- If plaster is used in place of the foot shoes.
- with adjustable Velcro tapes can also be applied to both feet.

S	34-37 cm	Color	Black
M	38-41 cm		
L	42-45 cm		
XL	46-49cm		
XXL	49-52 cm		

Ref-100.030



OrthoCy Hallux Valgus Support for Night

- Helps to keep the hallux (big toe) aligned in a desired position with the help of stay.
- Designed to use during night time.

Adjustable	Color	<input checked="" type="checkbox"/> Dark Blue
		<input type="checkbox"/> Grey

Ref-100.032



OrthoCy Hallux Valgus Support for Day

- Helps to keep the big toes aligned, designed to use during day time.
- Velcro and elastic bandage is made of neo prene fabric. Day (shoes, slippers ... in) designed to be used.

Adjustable	Color	<input checked="" type="checkbox"/> Black
		<input type="checkbox"/> Black
Color	S	34-37 cm
	M	38-41 cm
	L	42-45 cm
	XL	-
Color	XXL	-

Ref-105.246



OrthoCy Silicon Metatarsal Cushion

- Made of silicone material and the top surface is fabric-coated.
- Transversally, it is placed proximal to the metatarsal heads.
- Metatarsi reduces pain at the by relieving pressure.

Color	S	35-37 cm
	M	38-40 cm
	L	41-43 cm
	XL	-
Color		<input checked="" type="checkbox"/> Darkblue

Ref-105.247



OrthoCy Silicone Heel Cushion with Fabric

- It has relaxing features when you have in Epin calcaneus and achilles tendon.
- Heel Cushion With Silicone Fabric
- It has relaxing features when you have in Epin calcaneus and achilles tendon.

Color	S	35-37 cm
	M	38-40 cm
	L	41-43 cm
	XL	-
Color		<input checked="" type="checkbox"/> Darkblue

Ref-105.248



OrthoCy Internal-External Wedge

- Where inward or outward push the soles; shoes placed into the medial or lateral wedge surface (elevation) provides.
- Middle degree of hardness, is made from 100% silicon

Color	S	35-37 cm
	M	38-40 cm
	L	41-43 cm
	XL	-
Color		<input checked="" type="checkbox"/> Darkblue

Ref-105.249



OrthoCy Silicone Insoles with Fabric

- Provides anatomical support to the sole and enables even distribution of pressure on the sole.
- Has a relaxing and comforting effect.
- Middle degree of hardness, is made from 100% silicon.

Color	S	35-37 cm
	M	38-40 cm
	L	41-43 cm
	XL	43-45 cm
Color		<input checked="" type="checkbox"/> Darkblue

Ref-105.250



OrthoCy Silicone Insoles

- Provides anatomical support to the sole and enables even distribution of pressure on the sole.
- Has a relaxing and comforting effect.
- Middle degree of hardness, is made from 100% silicon

	S	M	L	XL	XXL
	35-37cm	38-40cm	41-43cm	43-45cm	45-47cm

Ref-105.251




OrthoCy Silicone Insoles 3 / 4

- Provides anatomical support to the sole and enables even distribution of pressure on the sole.
- Has a relaxing and comforting effect.
- Middle degree of hardness, is made from 100% silicon.
- 3/4 in shoes be in size and offers comfortable use.

	S	M	L	XL	XXL
	35-37cm	38-40cm	41-43cm	43-45cm	45-47cm

Ref-105.252



OrthoCy Silicone Heel Cushion

- Provides a relaxing effect in case of epin cal canus and achilles tendon irritations.
- Middle degree of hardness, is made from silicon

	S	M	L	XL
	35 -37	38 -40	41-43	44-46

Ref-105.253



OrthoCy Silicone Arch Support with Fabric

- Provides excellent support on the longitudinal arch in flat foot conditions.
- Low degree of hardness is made of silicone.
- Contact with the base portion is covered with soft fabric.
- Shoe has a bonding surface.

	S	M	L	XL
	35 -37	38 -40	41-43	-

Ref-105.254



OrthoCy Silicone Toe Separator

- Helps proper positioning of toes in the shoes and prevents overlapping.
- Low degree of hardness is made of 100% silicone.

	S	M	L	XL
	35 -37	38 -40	41-43	-

Ref-105.255



OrthoCy Metatarsal Pad

- Provides support for metatarsal curve in bio mechanical problems of metatarsal bones.
- Neoprene fabric with silicone pads built into the cozy and comfortable to use

	Adjustable
	Adjustable

Ref-105.260



OrthoCy Silicone Metatarsal Pad with Fabric

- Wedge is used in between the toes and supports metatarsal area. It uses to relaxing sensitive of metatarsal head and metatarsalgia sick.
- Low degree of hardness is made of 100% silicone.
- Wedge is used in between the toes and supports metatarsal area.
- It uses to relaxing sensitive of metatarsal head and metatarsalgia sick.

	Adjustable
	Adjustable

Ref-105.261



OrthoCy Silicone Metatarsal Pad with Ring

- Metatarsal is relaxed when the ring passed to second toe.
- It uses to relaxing sensitive of metatarsal head and metatarsalgia sick.

	Adjustable
	Adjustable

Ref-105.262



OrthoCy Silicone Bunion Protector

- The bunion protects to shoes pressure and help to relaxing callused area.
- It used in Hallux valgus treatment.



Adjustable

Ref-105.263



OrthoCy Silicone Cushion Hammer Toe

- It protect the thumb when you have hammer finger problem.
- It gives relaxing and prevent curling.



Adjustable

Ref-105.264



OrthoCy Finger Sleeve

- When you have toe friction, it makes callus.
- This products reduce to pressure and pain on callus.



	S	M	L	XL
	35 -37	38 -40	41-43	

Ref-105.267



OrthoCy Silicone Bunion Protector with Roller

- Spacer aligns the big and second toes to help relieve pain and stress on the bunion joint
- Medical-grade silicone gel.
- Comfortable Loop Keeps Spreader in Place
- One size fits most and easy to use
- Helps to prevent callous and corns from forming



Adjustable

Ref-106.243



OrthoCy Insoles (Heel Open)

- Thanks to the hole in the epin section, heal ing helps in heel area complaints.



35.....44

Ref-106.244



OrthoCy Insoles (Big Toe)

- Many nerves pass through the foot and these nerves are responsible for giving their feet feelings.
- The nerve passing between the fingers 3. and 4. can be under pressure due to the anatomical structure.



35.....44

Ref-106.245



OrthoCy Insoles (Arch Wedge)

- Provides arch support and ease of use on feet with heel zone.
- It is especially used in narrow or wide shoe types in women's heeled shoes.



35.....44

Ref-106.246



OrthoCy Insoles 3 / 4 (Hard Foam)

- Provides anatomical support to the base of the foot that the pressure on the base is balanced.
- The base is 3/4 in size, allowing comfortable use within the shoe.
- Made from Eva (medium hard) material.



35.....44

Ref-106.247



OrthoCy Insoles Arc Reinforced (Leather)

- The Latex cushion in the middle of the foot helps balance the foot load.
- It gives comfort to the foot.
- Thanks to the extra latex support on the heel, the pressure on the heel is reduced.



35.....44

Ref-106.248



OrthoCy Insoles Arch Reinforced (Eva Foam)

- Support to the anatomic sole of the foot base and provides balanced adjustment of the pressure in the foot base.
- Made of Eva (medium hardness) material.
- The foot base is designed to fill all the spaces.



35.....44

Ref-106.256



OrthoCy Insoles Flexible S / 4

- Provides excellent support on the longitudinal arch in flat foot conditions.
- Low degree of hardness is made of 100% silicone.
- Contact with the base portion is covered with soft fabric.
- Shoe has a bonding surface.



35.....44

Ref-106.257



OrthoCy Insoles (Steel)

- Provides firm support in the event of pes planus and pes cavus with its steel arch.
- Longitudinal and metatarsal support at the holes on the original leather covered steel mold.



35.....44

Ref-106.258



OrthoCy Insoles (Cork)

- Fits anatomically to the sole and provides ideal support for longitudinal arch and transverse arch.
- Longitudinal and metatarsal support at the holes on the original leather covered cork mold.



35.....44



Varicose Socks

The Combination of
comfort, technology,
design and functionality




Ref-110.220



VarCy Arm Compression (Shoulder Wrist)

- The treatment of congenital lymphoedema and acquired lymphoedema,
- For the arms which may become prone to swelling after removal of a cancerous breast and associated lymph nodes,
- Postoperative and post-traumatic oedema of the upper limbs,
- Incidental trauma such as burns and lacerations.



(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49


S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.221



VarCy Arm Compression (Shoulder Elbow)

- The treatment of congenital lymphoedema and acquired lymphoedema,
- For the arms which may become prone to swelling after removal of a cancerous breast and associated lymph nodes,
- Postoperative and post-traumatic oedema of the upper limbs,
- Incidental trauma such as burns and lacerations.



(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49


S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.431



VarCy Varicose Socks Lower Compressed (G-T Calf)

- Provides compression and blood circulation to the calf.
- Applying pressure of 14-20 mmHg.
- It offers excellent comfort and health in legs.
- It is made of 140 denier, very thin and knitting technique.



(SIZE TABLE)	cB	cD	cF	cG
1-Small	19-21	29-32	39-47	42-55
2-Medium	21-23	30-35	41-50	44-58
3-Large	23-25	32-38	45-54	48-62
4-X Large	25-27	33-42	46-57	50-67
5-XX Large	27-29	35-45	47-60	53-73


S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.432



VarCy Varicose Socks Lower Compressed (G-T Thigh)

- Provides compression and blood circulation to the calf.
- Applying pressure of 14-20 mmHg.
- It offers excellent comfort and health in legs.
- It is made of 140 denier, very thin and knitting technique.



(SIZE TABLE)	cB	cD	cF	cG
1-Small	19-21	29-32	39-47	42-55
2-Medium	21-23	30-35	41-50	44-58
3-Large	23-25	32-38	45-54	48-62
4-X Large	25-27	33-42	46-57	50-67
5-XX Large	27-29	35-45	47-60	53-73


S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.433



VarCy Varicose Socks Lower Compressed (G-T Panty)

- Provides compression and blood circulation to the leg.
- Applying pressure of 14-20 mmHg.
- It offers excellent comfort and health in legs.
- It is made of 140 denier, very thin and knitting technique.



(SIZE TABLE)	cB	cD	cF	cG
1-Small	19-21	29-32	39-47	42-55
2-Medium	21-23	30-35	41-50	44-58
3-Large	23-25	32-38	45-54	48-62
4-X Large	25-27	33-42	46-57	50-67
5-XX Large	27-29	35-45	47-60	53-73


S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.434



VarCy Varicose Socks Lower Compressed (G-T Pregnancy)

- Provides compression and blood circulation to the leg in case of varicose pains during pregnancy.
- Applying pressure of 14-20 mmHg.
- It offers excellent comfort and health in legs.
- It is made of 140 denier, very thin and knitting technique.



(SIZE TABLE)	cB	cD	cF	cG
1-Small	19-21	29-32	39-47	42-55
2-Medium	21-23	30-35	41-50	44-58
3-Large	23-25	32-38	45-54	48-62
4-X Large	25-27	33-42	46-57	50-67
5-XX Large	27-29	35-45	47-60	53-73

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.531



Nose Open - Off

VarCy Varicose Socks Medium Compressed (C-2 Calf)

- Used in the treatment of beginning and medium level varicose veins by increasing blood circulation.
- Applying pressure of 20-30 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 240 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.532



Nose Open - Off

VarCy Varicose Socks Medium Compressed (C-2 Thigh)

- Used in the treatment process of medium level varicose veins to increase blood circulation.
- Applying pressure of 20-30 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 240 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.533



Nose Open - Off

VarCy Varicose Socks Medium Compressed (C-2 Panty)

- Provides medium level compression and blood circulation to the leg.
- Used in the treatment process of medium level of varicose veins to increase blood circulation.
- Applying pressure of 20-30 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 240 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.534



Nose Open - Off

VarCy Varicose Socks Medium Compressed (C-2 Pregnancy)

- Provides compression and blood circulation to the leg in case of varicose pains during pregnancy.
- Used in the treatment process of medium level varicose veins to increase blood circulation.
- Applying pressure of 20-30 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 240 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.631



Nose Open - Off

VarCy Varicose Socks High Compressed (C-3 Calf)

- Provides high compression and blood circulation to the calf.
- Used in the treatment process of varicose veins.
- Applying pressure of 30-40 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 280 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.632



Nose Open - Off

VarCy Varicose Socks High Compressed (C-3 Thigh)

- Provides high compression and blood circulation to the calf and thigh.
- Used in the treatment process of varicose veins.
- Applying pressure of 30-40 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 280 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1-Small	16-19	13-15	21-27	23-30
2-Medium	18-21	15-17	24-30	26-33
3-Large	20-23	17-19	27-33	30-37
4-X Large	22-25	19-21	30-36	34-41
5-XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	Nude
1	2	3	4	5	

Ref-110.633



Nose Open - Off

VariCy Varicose Socks High Compressed (1-3 Party)

- Provides high compression and blood circulation to the leg.
- Used in the treatment process of varicose veins.
- Applying pressure of 30-40 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 280 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1 - Small	16-19	13-15	21-27	23-30
2 - Medium	18-21	15-17	24-30	26-33
3 - Large	20-23	17-19	27-33	30-37
4 - X Large	22-25	19-21	30-36	34-41
5 - XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	
1	2	3	4	5	Nude

Ref-110.634



VariCy Varicose Socks High Compressed (1-3 Pregnancy)

- Provides compression and blood circulation to the leg during pregnancy to prevent varicose pains .
- Applying pressure of 30-40 mmHg.
- It offers excellent comfort and health in legs.
- Available as open and close toe options.
- It is made of 280 denier, very thin and knitting technique.

(SIZE TABLE)	cA	cC	cE	cG
1 - Small	16-19	13-15	21-27	23-30
2 - Medium	18-21	15-17	24-30	26-33
3 - Large	20-23	17-19	27-33	30-37
4 - X Large	22-25	19-21	30-36	34-41
5 - XX Large	24-27	23-25	36-42	42-49

S	M	L	XL	XXL	
1	2	3	4	5	Nude

Ref-110.731



VariCy Anti-embolism Socks (Calf)

- Used in the recovery process of the patients to prevent thrombosis after the surgery and pre or after labour or bedridden patients to improve blood flow.
- The hole on the toe area helps to check the blood circulation in an easy way.
- It is made of 140 denier, very thin and knitting technique.

(SIZE TABLE)	cB	cD	cF	cG
1 - Small	19-21	29-32	39-47	42-55
2 - Medium	21-23	30-35	41-50	44-58
3 - Large	23-25	32-38	45-54	48-62
4 - X Large	25-27	33-42	46-57	50-67
5 - XX Large	27-29	35-45	47-60	53-73

S	M	L	XL	XXL	
1	2	3	4	5	White

Ref-110.732



VariCy Anti-embolism Socks (Thigh)

- Used in the recovery process of the patients to prevent thrombosis after the surgery and pre or after labour or bedridden patients to improve blood flow.
- The hole on the toe area helps to check the blood circulation in an easy way.
- The silicone band around the knee area prevents slipping.
- It is made of 140 denier, very thin and knitting technique.

(SIZE TABLE)	cB	cD	cF	cG
1 - Small	19-21	29-32	39-47	42-55
2 - Medium	21-23	30-35	41-50	44-58
3 - Large	23-25	32-38	45-54	48-62
4 - X Large	25-27	33-42	46-57	50-67
5 - XX Large	27-29	35-45	47-60	53-73

S	M	L	XL	XXL	
1	2	3	4	5	White

Ref-110.733





Right / Left


VariCy Anti-embolism Socks (Thigh - Waist Belt)

- Used in the recovery process of the patients to prevent thrombosis after the surgery and pre or after labour or bedridden patients to improve blood flow by providing calf compression.
- The hole on the toe area helps to check the blood circulation in an easy way.
- It is made of 140 denier, very thin and knitting technique.

(SIZE TABLE)	cB	cD	cF	cG
1 - Small	19-21	29-32	39-47	42-55
2 - Medium	21-23	30-35	41-50	44-58
3 - Large	23-25	32-38	45-54	48-62
4 - X Large	25-27	33-42	46-57	50-67
5 - XX Large	27-29	35-45	47-60	53-73

S	M	L	XL	XXL	
1	2	3	4	5	White





OneCy

One Size

The Combination of
comfort, technology,
design and functionality

Ref-120.103



OneCy Bandage Cervical Orthosis

- Provides limits cervical spine motion to varying degrees.
- Used in the treatment of a disc hernia cervical osteoarthritis and cervical muscle pains.
- Thanks to the bandage feature applicable to all neck structures.
- It is made of polyethylene material covered with sponge fabric laminated sponge.



■ Nude
■ Black

Adjustable Size

Ref-120.230



OneCy Corset Lumbosacral

- Provides support to lower back muscles in the treatment process of facet joint disorders.
- It is made of neoprene fabric, which has a velcro holding feature.

■ Black

Adjustable Size

Ref-120.420



OneCy Corset Abdominal

- Provides support to the abdominal area following surgeries and prevents sagging.
- It is made of fabric, which has a velcro holding feature.



■ Black

Adjustable Size

Ref-30.300



OneCy Arm Sling

- Stabilizes the arm in the desired position to immobilize the arm.
- It is made of fabric, which has a velcro holding feature.

■ Black

Adjustable Size



OneCy Neoprene Epicondylitis Bandage

- Used in the treatment of epicondylitis by decreasing the load on the tendons by applying compression.
- Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epicondyle heads.
- It is made of neoprene fabric, which has a velcro holding feature.

Ref-50.303



Adjustable Size

Black

Ref-120.210



OneCy Elbow Support

- Provides warmth and support for the elbow joint.
- Used in the treatment of tennis elbow syndrome and Epichondylitis.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

Ref-60.600



OneCy Neoprene Hand Strap

- Provides protective support for wrist during sports activities.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

Ref-60.603



OneCy Neoprene Wrist Bandage

- Provides protective support during sports activities.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

OneCy Wrist Support

- Provides immobilization and constrains to the movements of the wrist
- It is made of neoprene fabric, which has a velcro holding feature.

Ref-120.211



Adjustable Size

Black

OneCy Thumb Support

- Provides support to the wrist joint and stabilization in mild strains.
- Used in De quarvain syndrome and thumb ligament problems.
- It is made of neoprene fabric, which has a velcro holding feature.

Ref-120.213



Adjustable Size

Black

Ref-120.221



OneCy Wrist & Thumb Support

- Provides excellent immobilization to the thumb joint and gives an external support.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

OneCy Neoprene Knee Tendon Band

- Provides supports the patellar tendon with silicone pad.
- Used in the sports activities and treatment of Osgood –Schlatter's syndrome and patellar tendonitis.
- It is made of neoprene fabric, which has a velcro holding feature.

Ref-20.110



Adjustable Size

Black

Ref-120.200



OneCy Knee Support

- Used in treatment process of patellar femoral pains, osteoarthritis and meniscus injuries.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

Ref-120.205



OneCy Knee Support Patella Fixed

- Used in the treatment process of patellar femoral pains, osteoarthritis and meniscus injuries.
- Reinforced patella support gives high support to knee area.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

Ref-90.406



OneCy Neoprene Ankle Support & Bandage

- Provides warmth and support for ankles with mild pain.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black

Ref-120.220



OneCy Ankle Support

- Provides warmth and support for the ankle joints
- Provides stabilization in mild strains.
- It is made of neoprene fabric, which has a velcro holding feature.

Adjustable Size

Black



KiddoCy

Kiddo Product

The Combination of
comfort, technology,
desing and functionality



Ref-130.370



KiddoCy Cervical Orthosis (Sponge)

- The UCBL used in childhood, with its hard-core structure and special design, distributes the whole body equally to the feet.
- Support to the heel bone in a fully orthopedic manner with a gauge and a suitable design.



One size only

COLOR

- White
- Darkblue

Ref-130.300



KiddoCy Arm Sling

- Stabilizes the arm in the desired position. (for children)
- Adjustable neck strap available.
- Made of laminated sponge with terry cloth inner surface and outer surface part of velcro holding fabric.



Adjustable

COLOR

- Darkblue

Ref-130.301



KiddoCy Arm Sling (Mesh Fabric)

- Stabilizes the arm in the desired position to immobilize the arm-shoulder. (for children)
- Made of mesh fabric.



S	16-22 cm
M	22-28 cm
L	28-34 cm
XL	-
XXL	-

COLOR

- White

Ref-130.302



KiddoCy Arm Sling (Printed Fabric)

- Stabilizes the arm in the desired position to immobilize the arm-shoulder.
- Made of gabardine fabric special printed for children.

	S	16-22 cm	Color	■ DarkBlue
	M	22-28 cm		
	L	28-34 cm		
	XL	-		
	XXL	-		

Ref-130.311



KiddoCy Clavicle Bandage

- It can replace a plaster to stabilize the clavicle bone in cases of clavicle fractures.
- The one-piece adjustable bandage is designed- for the treatment of clavicle fractures.
- Made of laminated sponge in thickness.

	S	40-50 cm	Color	■ DarkBlue
	M	50-60 cm		
	L	60-70 cm		
	XL	-		
	XXL	-		

Ref-130.312



KiddoCy Clavicle Bandage

- It can replace a plaster to stabilize the clavicle bone in cases of clavicle fractures.
- The one-piece adjustable bandage is designed- for the treatment of clavicle fractures.
- Made of laminated sponge in thickness.



	S	40-50 cm	Color	■ DarkBlue ■ Red
	M	50-60 cm		
	L	60-70 cm		
	XL	-		
	XXL	-		

Ref-130.380




Ref-130.379



KiddoCy Joint Hyperextension

- Stabilizes the hyperextension of vertebral colon on the sagittal plan compression fractures of the vertebra.
- Inhibits lateral flexion with aluminum frame.



One size only

Ref-130.320



KiddoCy Elbow Support

- Provides warmth and support for the elbow joint. (for children)
- Used in the treatment process of tennis elbow syndrome and epichondylitis.
- Made of dual laminated perforated neoprene which allows air through the channels.

- Red
- Blue


	S	16-18 cm
	M	18-20 cm
	L	20-22 cm
	XL	-
	XXL	-

Ref-130.326



KiddoCy Wrist Bandage

- Provides protective support for wrist during sports activities. (for children)
- Made of dual laminated perforated neoprene which allows air through the channels.



- Red
- Blue

Adjustable Size

Ref-130.327



KiddoCy Wrist Support

- Provides immobilization and constrains to the movements of the wrist. (for children)
- Made of dual laminated perforated neoprene which allows air through the channels.

- Red
- Blue

	S	9-11 cm
	M	11-13 cm
	L	13-15 cm
	XL	-
	XXL	-

Left - Right

Ref-130.328



KiddoCy Wrist & Thumb Support

- Enables immobilization of the wrist in conditions such as carpal tunnel syndrome and painful wrist.
- The aluminum rod inside constrains the movements of the wrist.
- Made of dual laminated perforated neoprene which allows air through the channels.

- Red
- Blue

	S	9-11 cm
	M	11-13 cm
	L	13-15 cm
	XL	-
	XXL	-

Left - Right

Ref-130.330



KiddoCy Knee Immobilizer

- Provides strong stabilization and extension to the knee joint.
- Used in the treatment process of post patellar, arthral, and collateral ligament operations.
- Made of sponge fabric lined with a terry cloth on its inner surface.

	S	25-28 cm
	M	28-31 cm
	L	31-34 cm
	XL	-
	XXL	-

Ref-130.350



KiddoCy Knee Support

- Provides protection and warmth.
- Functions well with mild traumas and inflammation.
- Made of dual laminated perforated neoprene which allows air through the channels.

	S	25-28 cm
	M	28-31 cm
	L	31-34 cm
	XL	-
	XXL	-

Ref-130.351



KiddoCy Knee Support With Flexible Stays

- Supports collateral ligament with the flexible stays on the sides.
- Used in especially ligament injuries.
- Flexible steel springs on the both sides provide convenient support to the knee in case of collateral traumas.
- Made of dual laminated perforated neoprene which allows air through the channels.

	S	25-28 cm
	M	28-31 cm
	L	31-34 cm
	XL	-
	XXL	-

Ref-130.360



KiddoCy Ankle Support

- Provides warmth and support for the ankle joints.
- Provides stabilization in mild strains.
- Made of colored neopren fabric with thickness. (for children)

	S	18-20 cm
	M	20-22 cm
	L	22-24 cm
	XL	-
	XXL	-

Ref-130.361



KiddoCy Ankle Support (S-Bandage)

- Provides warmth and support to the ankle in case of moderate pain.
- Made of colored neopren fabric with thickness.

	S	18-20 cm
	M	20-22 cm
	L	22-24 cm
	XL	-
	XXL	-

Ref-130.381



KiddoCy Ankle Foot Orthosis (AFO)

- For children, helps to keep the ankle in functional position of 90°.
- Made of polyurethane in 3 mm thickness.
- Provides comfort and cleanliness thanks to removable terrycloth pads.



	S	21..24 A.No
	M	25..28 A.No
	L	28..33 A.No
	XL	-
	XXL	-



Ref-133.001

KiddoCy Dennis Brown Orthosis

- Dennis brown boots are specially produced for children used with pectus after hip release surgery.



Ref-133.002

KiddoCy Dennis Brown Orthosis with Straps


- Dennis brown boots specially produced for children are used with pectus after hip release surgery.
- Due to the special hole in the heel area, the fit between the foot and the boots can be easily observed.
- It is made from authentic leather material, does not irritate the feet.



Ref-130.245

KiddoCy UCBL Insoles


- The UCBL used in childhood, with its hard-core structure and special design, distributes the whole body equally to the feet.
- Support to the heel bone in a fully orthopedic manner with a gauge and a suitable design.



Ref-132.246

KiddoCy Insoles 3 / 4 (Hard Foam)

- Used in orthopedic foot problems (Pescavus, Pes planus, Plantar fasciitis).
- More active lifestyle; children, young and middle age groups and heavily overweight people are preferred to use 3/4 insoles.



Ref-130.259

KiddoCy Insoles Clark with Arch Reinforcement (Eva Foam)

- Provides support for longitudinal arch and transfer arch while anatomically conforming to the base of the foot.



Ref-130.422



KiddoCy Corset Hernia Support

- It pressure abdominal area and prevents the progression of hernia
- Thanks to special hernia of pad, when you have abdominal hernia.
- It made from elastic woven fabric with using cotton yarn.

S	03 Month
M	03-06 Month
L	06-09 Month
XL	-
XXL	-

Ref-132.003



KiddoCy GNC Hip Orthosis (Plastic)

- Used to maintain abduction and articulation of the head of the femur with the acetabulum in a baby born with dislocated hips.



S	03 Month
M	03 – 06 Month
L	06 – 09 Month
XL	-
XXL	-

Ref-132.001



KiddoCy GNC Hip Orthosis (Yveska)

- It pressure abdominal area and prevents the progression of hernia thanks to special hernia of pad, when you have abdominal hernia.
- It made from elastic woven fabric with using cotton yarn.




S	03 Month
M	03 - 06 Month
L	06 - 09 Month
XL	09 - 12 Month
XXL	12 - 15 Month

Ref-130.390



KiddoCy GNC Hip Orthosis (Pavlik Bandage)

- A Pavlik harness is used in the treatment of hip dysplasia in infants.
- The Pavlik harness is fit to the baby and holds his or her hips in proper position.
- By flexing up the legs, and allowing the knees to fall outwards, the hips are held in proper position.

	S	03 Month
	M	03 - 06 Month
	L	06 - 09 Month
	XL	-
	XXL	-



lifeCy

Life Supports

The Combination of
comfort, technology,
design and functionality



Ref-142.101

Viscy Pillow Back

- Provides comfort and Support to the lumbar area.
- Keeps your spine in a straight line.



Ref-142.110



Viscy Pillow Cushion Set

- Provides comfort and relaxing effect in case of coccyx traumas.



Ref-142.120

Viscy Pillow Neck

- Provides full support to the shoulder and neck muscles by providing ideal sleeping position.
- It is easy to cleaning with zippered cover.
- Made from water-based polyurethane material.



Ref-142.130

Viscy Pillow Reflex

- Provides comfortable sleep position in case of reflux problems.



Ref-142.125

Viscy Pillow Travel

- Provides comfort and support to the neck.



Ref-171.180

Mattress and Pillow Protectors (50*70 cm)

- It is suitable for longterm staying patients, bedwetting problem of children or adults, and people how have sweltering problem at night.
- It has air circulation features due to micro pore.
- No water exist, It prevents formation of bacteria, It is durable for spots and malodours.
- Upper layer is laminated by soft terry fabric and bottom layer is laminated by protective PU membrane fabric.
- It covers bed exactly, spills and non creasing due to elastic bandage in its corners.

Size 50X70 cm



Ref-171.190

Ref-171.200

Mattress and Pillow Protectors (90*190 cm)- (160*200 cm)

- It is suitable for longterm staying patients, bedwetting problem of children or adults, and people how have sweltering problem at night.
- It has air circulation features due to micro pore.
- No water exist, It prevents formation of bacteria, It is durable for spots and malodours.
- Upper layer is laminated by soft terry fabric and bottom layer is laminated by protective PU membrane fabric.
- It covers bed exactly, spills and non creasing due to elastic bandage in its corners.

Size 90 X 190 cm

Size 160 X 200 cm

Ref-105.230



Silicone Prosthesis Breast

- Silicone breast prosthesis is a real prosthesis that can follow all movements of your body and can be applied to both sides thanks to its symmetric structure.
- It is used after mastectomy operations.

Size
65-70-75-80-85-90-95-100-105-110-115-120-125-130

Ref-105.231



Bra for Breast Prosthesis

- The prosthesis is easily applied thanks to the pockets on both sides.
- The breast prosthesis is placed in these pockets and does not fall or slip.

Size
65-70-75-80-85-90-95-100-105-110-115-120-125-130

Ref-170.106

Ref-170.108

Ref-170.110

Ref-170.112

Ref-170.115

Ref-170.120



Elastic Bandage (8 cm) - (8 cm) - (10 cm) - (12 cm) - (15 cm) - (20 cm)

- In case of injuries in the body, in case of fracture or dislocation, protect the injured parts of the body by fixing the bandage.
- It is recommended to use the injured area or the wound in case of injuries because it protects against bacterial infection and bacterial infection.

Ref-154.027

Wheelchair Standard



Ref-154.030

Wheelchair (Functional)



Ref-154.101

Walker Aluminium



Ref-154.105

Walking Stick Aluminium



Ref-154.102

Crutches Aluminium (Pair)



Ref-154.103

Canadian Aluminium (Pair)



Ref-154.104

Triped Aluminium



Ref-154.106

Chaise Portable



 Morsa
Cyberg





Morsa Medikal Ortopedi ve Mak. San. Tic. Ltd. Şti.
İkitelli Organize Sanayi Bölgesi Giyim Sanayi Sitesi 4.Ada A Blok
No:602 Başakşehir / İstanbul / Türkiye
Tel: +90 212 671 11 39 Fax: +90 212 671 11 40
E-mail: info@morsa.com.tr - www.morsa.com.tr

 /MorsaCyberg
 @MorsaCyberg
 morsacyberg_ortopedi


Discover
the potential